

Tribeca

www.tribecafsp.com



Tribeca

HOSPITAL EQUIPMENTS



Emergency Trolley / Crash Cart


This product which is designed to keep all the medicines and medical devices for emergency code interventions, has the ergonomoy to allow more than one doctor and nurses for a fast and safe treatment at the same time around the trolley.



This product **colour alternatives, dimensions and detailed pictures** can be found at **Page 10 & 11.**



- 1 Serum Hanger**
Vertical Adjustable.
- 2 Defibrillator and Monitor Stand**
360 degree Rotatable.
- 3 Mechanical Breaking Lock**
Serial Number Plastic Seal
- 4 Nurse Writing Table**
Can be hidden when not in use.
- 5 Drawers**
3xh100mm short, 2xh150mm long drawer. Seperative trays inside the drawers can allow easy usage and optional medicine sorting..
- 6 Rack**
Mounting/Fitting Dock for dustbin, electrical plug and vacum pump.
- 7 Protective Bumper**
Your walls and trolley are safe by the bumper which is compatible to each corner.
- 8 Casters**
Thermo plactic plated wheels by 125mm caliber can allow more silent drive and easy rotation. Every wheels can be locked.tekerlek kilitlenebilir.



Tribeca

www.tribecafsp.com

9 Resuscitation Board

With 54cmx39cm dimension, lazer cut acrylic board.

10 Oxygen Tube Slot

110mm caliber, 3kg. tube with double strap connection, body docking station.

11 Probe Dock

Body docking station

12 File Shelves

Two-Piece Shelves

Address

ÜÇGEN ENDÜSTRİYEL MUTFAK ve OTEL EKİPMANLARI GIDA TURİZM SAN. TİC. LTD. ŞTİ.
ÜÇGEN INDUSTRIAL KITCHEN and OTEL EQUIPMENTS, FOOD, TOURISM, INDUSTRY and TRADE Ltd. Co.

Head Office: (For cargo and Warehouse): Akcaburgaz Mah. 4. Bolge 7. Cad. 124. Sok No:1 Kirac-Esenyurt – İstanbul 34522 TURKEY **Phone:** +90 212 886 71 91 - 90 **Fax:** +90 212 886 71 94

Branch Office: Aspendos Bulvari No: 204/4 07200 Antalya / TURKEY **Phone & Fax:** +0 242 311 94 34

Tribeca

www.tribecafsp.com