



Therapy, Medication and Anesthesia Trolley

Drawers with tray and divided separate system allows for high capacity storage and optional arrangement possibility. Both lightweight and ergonomic design and as well esthetic appearance designed to fulfill all requirement of the professional corporations.



This product **colour alternatives, dimensions and detailed pictures** can be found at **Page 10 & 11.**



- 1 Serum Hanger**
Vertical Adjustable
- 2 Tilting Box Cover**
10,12,18 units by 3 different size option
- 3 Key With Lock Mechanism**
Serial Number Plastic Seal.Optional Key Lock
- 4 Nurse Writing Table**
Can be hidden when not in use.
- 5 Drawers**
3x100mm short, 2x150mm long drawer. Seperate trays inside the drawers can allow easy usage and optional medicine sorting.7 colour alternatives anf maximum 6 drawer usage.
- 6 Waste Bin**
Body mounted, tilting cover waste bin can be taken out and easily cleaned in case of need.
- 7 Protective Bumper**
Your walls and trolley are safe by the bumper which is compatible to each corner.

8 Casters

Thermo plastic plated wheels by 125mm caliber can allow more silent drive and easy rotation. Every wheels can be locked.

File Shelves*

Two-Piece Shelves

Oxygen Tube Slot*

110mm caliber, 3kg. tube with double strap connection, body docking station.

Probe Dock*

Body docking station

* The visuals of these features can be seen at **page 5**.

Tribeca

www.tribecafsp.com

