



MY-PRESSURE 2.0 is an electronic digital device for the measurement of blood pressure and heart pulse. Provides easy and accurate reading of systolic and diastolic pressure, heartbeat and possible arrhythmia. Fully automatic operation along with the large monitor makes the device easy to use and simple to read. The software for measurement of MYPRESSURE 2.0 is clinically validated and in compliance with the ESH protocol (European Society of Hypertension).

OPERATING FEATURES

- Fuzzy logic technology
- Irregular heartbeat indicator
- Indicator with colour scale in compliance with 1999 WHO-ISH guidelines
- Suitable for use with power adaptor (not included)
- 2 areas with 90 memories each (double patient)
- Average of the latest 3 measurements
- Date and Time function
- Automatic stand-by mode

STANDARD CONSUMABLES INCLUDED

- 1 x Cuff 22-32 cm circumference
- 1 x Carry Bag
- 4 x Alkaline Batteries (AA-type, 1.5V)

OPTIONAL ACCESSORIES

- 230V-50Hz Power Adaptor SP 0210
- Cuff 32-45 cm circumference SP 0197/01

TECHNICAL FEATURES

Size	cm 10,5 x 12 x 5,5
Weight	310 g
Measuring Method	Oscillometric
Inflation / Deflation Method	Automatic
Pressure Measurement Range	40-260 mmHg
Cuff Pressure Measurement Range	0-300 mmHg
Pulse/min (Pulse Rate)	40-160 beats/min
Pulse Rate Accuracy	+/- 5% Heart Rate Accuracy
Static Pressure Accuracy	+/- 3 mmHg
Power Source	4 AA-type Alkaline Batteries 1.5V (AC/DC Adaptor not included)
Years of Warranty	2
Units per Shipping Carton	10

CODE

RE 503000 MY-PRESSURE 2.0

CE 0123